

# Adjustable Loop BUTTON



## BUTTON LIGAMENT FIXATION SYSTEMS

### ADJUSTABLE BUTTON

#### ACL fixation systems

##### NEOSYS

offers a wide range of fixation devices dedicated to minimal invasive anterior cruciate ligament reconstruction surgeries (button with pre-measured loop size or adjustable version), repair of acromion

clavicular joint dislocation, fracture of the lateral quarter of the clavicle and Latarjet procedures.



SHOULDER



KNEE



#### FEMORAL FIXATION BUTTON CL

Pre-measured loop size button is made of :

- 1 titanium button (dimensions 12x4 mm)
- Polyester PET braided pre-measured loop
- 2 pre mounted UHMWPE sutures (green and white)

##### Indications

Tendon or ligament fixation in reconstruction surgeries especially for ACL repair.

##### Loop sizes

6 sizes are available and meet requirements of conventional ligamentoplasty surgeries (15, 20, 25, 30, 35, 40mm)





## ADJUSTABLE BUTTON LOOP (KNEE)

Adjustable button is made of :

- 1 titanium button (dimensions 12x4 mm)
- 1 UHMWPE braided double loop
- 2 pre mounted UHMWPE pulling wire (white and white/black)

### Indications

Tendon or ligament fixation in reconstruction surgeries and especially for ACL ligamentoplasty.



## ADJUSTABLE BUTTON LOOP (SHOULDER)

System is made of two devices :

- 2 adjustable button loop external (titanium button, UHMWPE double loop and pulling wires )

### Indications

Repair of acromium clavicular joint dislocation, lateral quarter of the clavicle fracture, Latarjet procedures.



**BUTTON LIGAMENT FIXATION  
SYSTEM**

**STORAGE**

**NEOSYS BUTTON** devices must be stored in original sealed packaging, in dry and clean environment, at ambient temperature.  
For any further information, please read the instructions for use.

**PRODUCTS**

**BUTTON** devices are delivered ready to use, sterile and packed in double sealed pouch.  
Sterilization method is EO according to the Norm ISO 11135  
Single use. **Do not re-use or re-process.**

Instruments for ACL/PCL ligamentoplasty, acromio-clavicular dislocation, Latarjet are **available upon request**.  
Do not hesitate to contact us for more details.

# NeoSys

## Medical Solutions

**[WWW.NEOSYSMED.COM](http://WWW.NEOSYSMED.COM)**

CE 0653  
CLASS IIb

