



## ARTOFIX Meniscal Repair System

Minimally Invasive | Simplicity | Efficiency

Meniscal Repair System


Meniscal Repair System consist of two PEEK implants and #2-0 UHMWPE Wilsuture, which provide a strong & simple minimally invasive solution for meniscal repair



All-inside Repair System

- ◇ **Minimally Invasive**  
Smaller implants and WilSuture Ultra #2-0 minimize disruption to the meniscus
- ◇ **Convenience**  
360° active deploy implants at any hand position.  
Auditory confirmation assures that implants are deployed
- ◇ **Detail**  
Adjustable depth limiter eliminates the risk of early deployment and damaging nearby neurovascular structures
- ◇ **High Efficiency**  
Reinforced needle rod penetrates with ease

## Meniscal Repair System



## All-inside Repair System

- ◇ **Convenience**
  - T1 passive deployed, T2 active deployed
  - Pre-tied sliding knot design eliminates the need of intra-articular knot tying
- ◇ **Minimal Invasive**
  - Small PEEK implants and WiSuture Ultra #2-0 ensure a strong and reliable meniscal repair
- ◇ **Versatility**
  - Tips are available in curved 17°, reverse 5° and straight, which helps access to different zone of meniscus.
- ◇ **High Efficiency**
  - Reinforced needle guarantees secure deployment



