


ORTHOSYN EGELOCK



REF	ÜRÜN ADI
16.150.024	Egelock ACL Femoral Tespit 15mm
16.150.025	Egelock ACL Femoral Tespit 20mm
16.150.026	Egelock ACL Femoral Tespit 25mm
16.150.027	Egelock ACL Femoral Tespit 30mm
16.150.028	Egelock ACL Femoral Tespit 35mm
16.150.029	Egelock ACL Femoral Tespit 40mm

Ön çapraz bağ cerrahisinde, femoral tünel hamstring greftin kortekse asılarak tespiti için tasarlanmış, etkinliği ve güvenilirliği klinik olarak kanıtlanmış bu düğme implant, ortopedi cerrahları tarafından en çok tercih edilen üründür.

Implant titanyumdan imal edilmiştir. Greftin asılacağı halka polyester örgüdür.

Çekme süturları işlem sırasında karışıklılığa yol açmamak için farklı renklerde dir.



ORTHOSYN EGELOCK

This button implant, which is clinically proven in effectiveness and reliability, is designed by orthopedic surgeons for anterior cruciate ligament surgery, femoral tunnel hamstring graft cortex hanging detection.

The implant is made of titanium. The loop is made of polyester woven .

Pulling sutures are different colors to avoid confusion during processing.

REF	PRODUCT DESCRIPTION
16.150.024	Egelock ACL button 15mm
16.150.025	Egelock ACL button 20mm
16.150.026	Egelock ACL button 25mm
16.150.027	Egelock ACL button 30mm
16.150.028	Egelock ACL button 35mm
16.150.029	Egelock ACL button 40mm



**NeoSys**
Cerrahi Çözümler
Address : 1322. Cadde 14/6 Öveçler, Ankara / Türkiye
Telephone : +90 312 479 0906 Fax : +90 850 226 0502
tbolat@neosysmed.com / export@neosysmed.com
<http://www.neosysmed.com>