


ARTHROSCOPY

Product features

- Endoscope adopts imported stainless steel material
- Endoscope adopts German optic glass, light fiber and light cone
- Adopt the patented technology of rod optical lens, clear field of view
- With direction index, sapphire lens cover never abrasion
- Stainless steel valve, easy for maintenance, avoid damaging
- Handle wheel of working element has ups and downs indicator
- Low temperature plasma sterilization can be used



Endoscope



CN020101

0° Endoscope
φ4mmx175mm



CN020102

30° Endoscope
φ4mmx175mm



CN020103

70° Endoscope
φ4mmx175mm

1. Connecting to STORZ light cable
Connecting to OLYMPUS light cable
2. Connecting to WOLF light cable
3. Connecting to ACMI light cable

NeoSys

www.neosysmed.com

ARTHROSCOPY

Sheath



CG070101

0° Sheath two-stopcock(Rotatable)
φ5.5mmx140mm



CG070201

30° Sheath two-stopcock(Rotatable)
φ5.5mmx140mm



CG070301

70° Sheath two-stopcock(Rotatable)
φ5.5mmx140mm

Trocar Cannula (Trocar)



CG060101

Prism puncture needed
φ4.3mmx184mm



CG060102

NeoSys










www.neosysmed.com


ARTHROSCOPY

- Multiple handle options to meet every Surgeon's preference and for optimized use of jaw
- Uniform and accurate force distribution
- Wide range of punches and jaw patterns
- Easy and precise use
- Cleaning through luer lock flush port on tube shaft instruments
- Maximum power transmission through a pinless joint guide
- Long durability due to the use of high-quality stainless steels and polymers
- Precise cutting pattern due to optimized punching gap dimensions in the jaw area



For Options of Punches based on the Application

	 Small Joint Arthroscopy	 Knee Arthroscopy	 Shoulder Arthroscopy	 Hip Arthroscopy
Working length	90mm	120mm/130mm	155mm/160mm	185mm
Classic handl 	●	●	●	●
Ergo handle 	n/a	●	n/a	●
Handle with positioning 			●	
Handle Slide-lock 	●	●	●	●
Cylinder handle handle 	n/a	n/a	●	n/a




NeoSys

www.neosysmed.com

ARTHROSCOPY

Oval Punch



Description	Bite Width	Tip Width	Ref#
straight	4.0 mm	5.80 mm	CY060521-01
upswept 15°	4.0 mm	5.80 mm	CY060521-02
curved left	4.0 mm	5.80 mm	CY060521-03
curved right	4.0 mm	5.80 mm	CY060521-0

Basket Punch Duckling



Description	Bite Width	Tip Width	Ref#
straight	3.2 mm	5.0 mm	CY061031-01
upswept 15°	3.2 mm	5.0 mm	CY061031-02
straight left	3.2 mm	5.0 mm	CY061031-03
straight right	3.2 mm	5.0 mm	CY061031-04
left curved, upward 15°	3.2 mm	5.0 mm	CY061031-05

Basket Punch Rotary



Description	Bite Width	Tip Width	Ref#
left90°	3.3 mm	7.9 mm	CY060331-01
right90°	3.3 mm	7.9 mm	CY060331-02


ARTHROSCOPY

Basket Punch Standard



Description	Bite Width	Tip Width	Ref#
straight	1.7 mm	4.1 mm	CY060821-01
upswept 15°	1.7 mm	4.1 mm	CY060821-02
left curved, upward 30°	1.7 mm	4.1 mm	CY060821-03
right curved, upward 30°	1.7 mm	4.1 mm	CY060821-04

Basket Punch



Description	Bite Width	Tip Width	Ref#
straight	2.35 mm	3.56 mm	CY080111-01
upswept 15°	2.35 mm	3.80 mm	CY080111-02
curved left	2.35 mm	3.56 mm	CY080111-03
curved right	2.35 mm	3.56 mm	CY080111-04

Reverse Basket Punch




Description	Bite Width	Tip Width	Ref#
left	2.3 mm	5.5 mm	CY060101-01
right	2.3 mm	5.5 mm	CY060102-01

NeoSys

www.neosysmed.com


ARTHROSCOPY

Suture Grasper




Tip Width	Working length	Ref#
4.0 mm	130mm/160mm/185mm	CY060205

Suture Grasper 1.0L



Tip Width	Tip Length	Working length	Ref#
5.2 mm	10.5 mm	130mm/160mm/185mm	CY030003

Suture Grasper 2.1D



Tip Width	Tip Length	Working length	Ref#
5.9 mm	6.7 mm	130mm/customized	CY030004-01


Suture Grasper 3.4L



Tip Width	Tip Length	Working length	Ref#
4.7 mm	5.3 mm	130mm/customized	CY030005-01


ARTHROSCOPY

Suture Grasper 11.3L




Tip Width	Tip Length	Working length	Ref#
5.3 mm	12 mm	130mm/customized	CY030002-01

Suture Grasper 1.0L



Tip Width	Working length	Ref#
4.1 mm	130mm/160mm/185mm	CY060203-01

Loose Body Forceps



Tip Width	Tip Length	Ref#
4.1mm	6.8 mm	CY060202-01


Cupped Grasper



Tip Width	Tip Length	Working length	Ref#
4.0 mm	7.3 mm	130mm/customized	CY060201-01


ARTHROSCOPY

Alligator Mouth Grasping Forceps



Tip Width	Working length	Ref#
4.0 mm	185mm/customized	CY030008-01

Suture Retrievers



Tip Width	Working length	Ref#
4.0 mm	90mm/160mm/185mm	CY060204-01



ARTHROSCOPY

Bird Beak



Description	Tip Width	Ref#
0°	3.4 mm	CY040004-01
Upswept10°	3.4 mm	CY040007-01
Upswept22°	3.4 mm	CY040002-01
Upswept30°	3.4 mm	CY040001-01
Upswept45°	3.4 mm	CY040003-01

Rhino Suture Passer Leftcurved



Description	Tip Width	Ref#
curved left	5.3 mm	CY040005


Rhino Suture Passer Rightcurved



Description	Tip Width	Ref#
curved right	5.3 mm	CY040006


ARTHROSCOPY

Suture Cutter



Bite Widht	Tip Width	Ref#
3.0 mm	5.0 mm	CY020001
2.7 mm	4.2 mm	CY020003

Suture Cutter



Bite Widht	Tip Width	Ref#
2.7 mm	5.50 mm	CY020002

Siliding Suture Cutter



Bite Widht	Tip Width	Ref#
4.0 mm	0.9 mm	CY020005-02



ARTHROSCOPY

Rasps



	Tip Widht	Working length	Ref#
Banana Knife	3.0 mm	105 mm	CG010101
Rosette Knife	3.0 mm	105 mm	CG010301
Bone File	3.0 mm	105 mm	CG010201

Right Angle Probe



	Tip Widht	Working length	Ref#
	3.0 mm	105 mm	CG040101
	5.0 mm	105 mm	CG040201

Curette



	Tip Widht	Working length	Ref#
Curette With Holes	3.0 mm	105 mm	CG050101
Curette	3.0 mm	105 mm	CG050201

ARTHROSCOPY

Knives



	Tip Widht	Working length	Ref#
Meniscus Knife	3.0 mm	105 mm	CG020101
Half-Round Knife	3.0 mm	105 mm	CG020201
Spade Knife	3.0 mm	105 mm	CG020301
Backward Knife	3.0 mm	105 mm	CG030101
Hook Knife	3.0 mm	105 mm	CG030201



ARTHROSCOPY

Knot Pusher



Tip Widht	Working length	Ref#
0.9 mm	165mm	CY020006-05
1.3 mm	165mm	CY020006-04
1.4 mm	160mm/180mm/185mm	CY020006-01-
2.0 mm	210mm	CY020006-06

Knot Pusher



Tip Widht	Working length	Ref#
0.3 mm	140mm/customized	CY020004


Knot Pusher



Tip Widht	Working length	Ref#
0.7 mm	145mm/customized	CY020007

ARTHROSCOPY


Knot Pusher



Description	Inside diameter	Tip Width	Ref#
right angle	1.0 mm	6.0 mm	CY050001
10° crescent	1.0 mm	6.0 mm	CY050004
curved left 20°	1.0 mm	6.0 mm	CY050005
curved right	1.0 mm	6.0 mm	CY050006
curved left 45°	1.0 mm	6.0 mm	CY050002
curved right 45°	1.0 mm	6.0 mm	CY050003
curved left 60°	1.0 mm	6.0 mm	CY050004
curved right 60°	1.0 mm	6.0 mm	CY020006-06
curved left 90°	1.0 mm	6.0 mm	CY020006-06

NeoSys

www.neosysmed.com



ARTHROSCOPY

- Especially designed to restore the rot. cuff and labral repair
- Use in combination with a single-use nitinol needle for
- Optimal suture (size #2) guidance in the operating room
- Slim line handle with optimized weight
- Optimized cleaning through luer lock flush port
- Easy puncturing of the tissue
- Available with and without finger ratchet
- Long finger handle for ergonomic safety activation of the jaws
- Trapdoor in top jaws for clamping the suture



Suture Needle 207 CY015001



Suture Needle 239 CY015002



Glenoid Labrum



Tip Width	Working length	Ref#
6.4 mm	170 mm	CY010002

Rotator Ruff



Tip Width	Working length	Ref#
6.6 mm	170 mm	CY010003