


# DISPOSABLE OATS SET

## DEVICE DESCRIPTION

- The Osteochondral Autograft Transfer System (OATS) includes the original OATS Set and the Single Use OATS Set.
- The OATS sets are size-specific cylindrical Harvesters with Collared Pins and various instruments to assist with graft insertion.
- The set includes a Pin, along with the optional Graft Driver and Core Extractor.
- It also includes a transparent Graft Delivery Tube for visualization of the graft during insertion.
- The Coring Reamers and Bone Graft Harvester sets are cylindrical bone core harvesters that include a collared pin.
- Suture Pass will simplify suture manipulation in arthroscopic procedures.

## INDICATIONS

- The OATS system is designed for Osteochondral Autograft Transplantation.
- The Coring Reamer System is designed to harvest a cylinder of cancellous bone while simultaneously creating the tibial tunnel.
- The harvested bone core can then be used to fill the patellar tendon harvest site, construct a bone-hamstring tendon-bone allograft, or enhance biological fixation of ACL grafts.
- The Bone Graft Harvester is designed to be used for various autograft bone harvesting procedures.
- Suture Pass will simplify suture manipulation in arthroscopic procedures.



## Order Informations

Ref Code	AR-1981-06S OATS SET 6.0 MM
Ref Code	AR-1981-0S OATS SET 8.0 MM
Ref Code	AR-1981-10S OATS SET 10.0 MM

NeoSys

[www.neosysmed.com](http://www.neosysmed.com)