

ORTHOSYN LIGAMENT STAPLE



REF	ÜRÜN ADI
16.150.030	Orthosyn Ligament Staple 8x20mm

Bağ, kemik veya tendon fiksasyonu için endikedir.

Staple çakıcı ve staple uyumu stabildir.

Maksimum kemik tutunumu için

güçlendirilmiş köşeler

Etkili anti-geri tepme tırtıkları.

Keskin eğik noktalar

Tüm anatomik veya protez bileşenlerin fiksasyonu için kullanılır.

Diz veya diğer eklemler için kullanılır.

Kobalt-krom alaşımı

İmplantın mikro

hareketinden dolayı ağrı riski azalır.

Steril



ORTHOSYN LIGAMENT STAPLE

Ligament, bone or tendon fixation
Full stability in the staple holder
Reinforced corners for stability.
Effective anti-backout serrations.

Sharp bevelled points
Used for the fixation of all anatomical or prosthetic components.
Used for knee or other joints
Cobalt-chromium alloy
Staples are fixed and stable in all directions.
Full stability in the staple holder ensuring immediate fixation stability.
Reduced risk of pain due to implant micromotion.
The staples are easy to remove because they provide a firm grip to the staple extractor.
Sterile

REF	PRODUCT DESCRIPTION
16.150.030	Orthosyn Ligament Staple 8x20mm



Address : 1322. Cadde 14/6 Öveçler, Ankara / Türkiye
Telephone : +90 312 479 0906 Fax : +90 850 226 0502
tbolat@neosysmed.com / export@neosysmed.com
<http://www.neosysmed.com>